

## Customized E2 Screening Kit (Cat. # J1200)

### Description

The Customized E2 Screening Kit (Cat# J1100) contains ubiquitin (Ub)-activating enzyme (Ube1), Ub, and 10 customer-selected Ub-conjugating enzymes (E2s) from our 29-listed E2s. It can be used for 1) identification of the E2 enzyme paired for a specific Ub ligase (E3) for a protein substrate, and 2) assessing/profiling E2 activity by monitoring the E2-Ub thioester conjugate formation under non-reducing condition.

### List of E2s and their reactivity (Select 10 E2s from the list below)

	E2 Name	Reactivity
1	6XHis-UbE2A	Ub
2	6XHis-UbE2B	Ub
3	6XHis-UbE2C	Ub
4	6XHis-UbE2D1	Ub
5	6XHis-UbE2D2	Ub
6	6XHis-UbE2D3	Ub
7	6XHis-UbE2D4	Ub
8	6XHis-UbE2E1	Ub
9	6XHis-UbE2E2	Ub
10	6XHis-UbE2E3	Ub
11	6XHis-UbE2F	Nedd8
12	6XHis-UbE2G1	Ub
13	6XHis-UbE2G2	Ub
14	6XHis-UbE2H	Ub
15	6XHis-UbE2I	SUMO

	E2 Name	Reactivity
16	6XHis-UbE2J1 <sub>(1-282)</sub>	Ub
17	GST-UbE2K	Ub
18	6XHis-UbE2L3	Ub
19	6XHis-UbE2L6	ISG15, Ub
20	6XHis-UbE2M	Nedd8
21	6XHis-UbE2N	Ub
22	6XHis-UbE2Q2	Ub
23	6XHis-UbE2R1	Ub
24	6XHis-UbE2R2	Ub
25	6XHis-UbE2S	Ub
26	6XHis-UbE2T	Ub
27	6xHis-Ubc13/UbE2V2	Ub
28	6XHis-UbE2W	Ub
29	6XHis-UbE2Z	Ub, FAT10
	UbE2Z works specifically with UBA6 E1.	

### Components

- 20X UBE1 (2µM) 50µl
- 10X E2 Ub-conjugating enzymes  
(Select 10 E2s from the list above, 20µM) 20µl each
- 10X Human Ubiquitin (500µM) 100µl
- 10X Ubiquitination Buffer 500µl
- 40 mM ATP 125µl

### Note



1. Reaction conditions should be optimized for specific assays. We recommend an initial testing condition as the following: a 20  $\mu$ l ubiquitination reaction contains 100nM E1, 2 $\mu$ M E2, 2 $\mu$ M E3, 2 $\mu$ M substrate, 50 $\mu$ M Ub, 2mM ATP, 1 $\mu$ l glycerol and 2 $\mu$ l 10X Ubiquitination Buffer. Substrate ubiquitination can be assayed by immunoblotting. According to this setup, the provided UbE1 and Ub are sufficient for 50 reactions; each E2 enzyme is enough for 10 reactions.
2. For monitoring the thioester conjugate formation between E2 and Ub to assess E2 activity, run non-reducing SDS-PAGE to preserve the thioester bond-linked conjugates.
3. 10X Ubiquitination Buffer: 200mM Tris, pH7.6 at 4°C, 500mM NaCl, 10mM  $\beta$ ME and 50mM MgCl<sub>2</sub>.
4. Store all components at -80°C upon receiving. Avoid multiple freeze-thaws.

