

## E2 Screening Kit II (Cat# J1300)

### Description

The E2 Screening Kit II (Cat# J1300) contains ubiquitin (Ub)-activating enzyme (Ube1), SUMO-activating enzyme (6xHis-SAE1/SAE2), Nedd8-activating enzyme (APPBP1/6xHis-UBA3), Ub, SUMO2, Nedd8, and 29 Ub-conjugating enzymes (E2s). It can be used for 1) identification of the E2 enzyme paired for a specific Ub ligase (E3), and 2) profiling E2 activity by monitoring the E2-Ub thioester conjugate formation under a non-reducing condition.

Comparing to the E2 Screening Kit (Cat# J1100), this kit contains additional SUMO or Nedd8-activating E1, SUMO2 and Nedd8, which allows to assess SUMOylation and neddylation.

### List of E2s and their reactivity

	E2 Name	Reactivity
1	6XHis-UbE2A	Ub
2	6XHis-UbE2B	Ub
3	6XHis-UbE2C	Ub
4	6XHis-UbE2D1	Ub
5	6XHis-UbE2D2	Ub
6	6XHis-UbE2D3	Ub
7	6XHis-UbE2D4	Ub
8	6XHis-UbE2E1	Ub
9	6XHis-UbE2E2	Ub
10	6XHis-UbE2E3	Ub
11	6XHis-UbE2F	Nedd8
12	6XHis-UbE2G1	Ub
13	6XHis-UbE2G2	Ub
14	6XHis-UbE2H	Ub
15	6XHis-UbE2I	SUMO

	E2 Name	Reactivity
16	6XHis-UbE2J1 <sub>(1-282)</sub>	Ub
17	GST-UbE2K	Ub
18	6XHis-UbE2L3	Ub
19	6XHis-UbE2L6	Ub, ISG15
20	6XHis-UbE2M	Nedd8
21	6XHis-UbE2N	Ub
22	6XHis-UbE2Q2	Ub
23	6XHis-UbE2R1	Ub
24	6XHis-UbE2R2	Ub
25	6XHis-UbE2S	Ub
26	6XHis-UbE2T	Ub
27	6xHis-Ubc13/UbE2V2	Ub
28	6XHis-UbE2W	Ub
29	6XHis-UbE2Z	Ub, FAT10
	UbE2Z works specifically with UBA6 E1.	

### Components

- 20X UBE1 (2µM) 100µl
- 20X 6xHis-SAE1/SAE2 (2µM, SUMO E1) 25µl
- 20X APPBP1/6xHis-UBA3 (2µM, Nedd8 E1) 25µl
- 10X E2 Ub-conjugating enzymes  
(29 E2s, see table above, 20µM each) 20µl each
- 10X Human Ubiquitin (500µM) 200µl
- 10X SUMO2 (500µM) 25µl
- 10X 6XHis-Nedd8 (150 µM) 50µl
- 10X Ubiquitination Buffer 1ml
- 20X ATP (40mM) 250µl

## Notes

1. Reaction conditions should be optimized for specific assays. We recommend an initial testing condition as the following for detecting protein ubiquitination: a 20 $\mu$ l ubiquitination reaction contains 100nM E1, 2 $\mu$ M E2, 2 $\mu$ M E3, 2 $\mu$ M substrate protein, 50 $\mu$ M Ub/50 $\mu$ M SUMO2/15 $\mu$ M Nedd8, 2mM ATP, 10% glycerol and 1X Ubiquitination Buffer (see table below and follow the addition order to mix components, tap the tube to mix well once water, ubiquitination buffer and glycerol are added). A negative control reaction can be included by omitting ATP in reaction mixtures. Incubate all reactions in a 37 °C water bath for 2 hours. Stop reactions by mixing with reducing SDS sample buffer. Protein ubiquitination can be detected by immunoblotting or Coomassie staining after SDS-PAGE separation. A similar reaction can be used to detect E3 autoubiquitination without a protein substrate.

Components	Reaction Concentration	Volume to be added for a 20 $\mu$ l reaction
Purified water	Not provided	6 $\mu$ l
10X Ubiquitination Buffer	1X	2 $\mu$ l
Glycerol	10% (not provided)	2 $\mu$ l
10X Human Ub, SUMO2 or Nedd8	50 $\mu$ M for Ub/SUMO2 15 $\mu$ M for Nedd8	2 $\mu$ l
10X Protein Substrate	2 $\mu$ M (not provided)	2 $\mu$ l
10X E2 Ub-conjug. Enzyme	2 $\mu$ M	2 $\mu$ l
10X E3 Ub Ligase	2 $\mu$ M (not provided)	2 $\mu$ l
20X UbE1	100nM	1 $\mu$ l
20X ATP	2mM	1 $\mu$ l

2. For monitoring the thioester conjugate formation between an E2 and Ub/SUMO2/Nedd8, a 20  $\mu$ l Ub charging reaction contains 100nM E1, 2 $\mu$ M E2, 50 $\mu$ M Ub/SUMO2 or 15 $\mu$ M Nedd8, 2mM ATP, 10% glycerol and 1X Ubiquitination Buffer. Non-reducing sample buffer should be used to stop the reaction and preserve the thioester bond-linked conjugates. Samples can be heated at 90 °C for 5 minutes. E2-Ub conjugates can be detected by immunoblotting or Coomassie staining after SDS-PAGE separation.
3. 10X Ubiquitination Buffer: 200mM Tris, pH7.6 at 4°C, 500mM NaCl, 10mM  $\beta$ ME and 50mM MgCl<sub>2</sub>. This buffer is suitable for ubiquitination, SUMOylation and neddylation reactions.
4. Store all components at -80°C upon receiving. Avoid multiple freeze-thaws.