

# Nickel XPure Agarose Resin

Cat. # P3010-5, P3010-25, P3010-100

**Size:** 5 ml (Cat.# P3010-5), 25 ml (Cat.# P3010-25), 100 ml (Cat.# P3010-100),

**Bead (Geometry, size):** Spherical, fine: ~ 20-50  $\mu\text{m}$

**Cross-Linked:** Yes

**Ligand:** Iminodiacetic acid (IDA)

**Agarose %:** 6% Agarose

**Binding Capacity:** > 26 mg/ml gel (Binding capacity will vary for each target protein.)


**Volume %:** 50% (v/v) aqueous suspension containing 20% Ethanol

**Application:** Batch, Gravity, MPLC and FPLC

**Introduction:** Nickel chelate resins recognize two exposed histidines for an appreciable retention of a protein with a 4X or 6XHis tag. Nickel XPure Agrose Resin provides the advantages of the Nickel cation with the high flow rates of the highly cross-linked agarose beads. Except for batch purification, this product is specially suitable for large scale His-tagged protein purifications using FPLC.

The XPure series of resins are made for large scale and fine purification using a high-performance liquid chromatography (HPLC) system.

**Image(s):**



Ube2I was purified with Nickle XPure Agarose Resin using a gravity column

Lane 1: Molecular weight marker  
Lane 2: Cell pellet after sonication  
Lane 3: Supernatant after sonication  
Lane 4: Flow through after binding  
Lane 5: Wash  
Lane 6: Elution 1  
Lane 7: Elution 2  
Lane 8: Elution 3  
Lane 9: Elution 4  
Lane 10: Elution 5

**Storage Temperature:** 2-8 °C

