


GST

Cat. # I1100

Also Known as: Glutathione S-transferases
NCBI Reference: N/A
MW (no tag): 28.2 kDa
Species: Human
Source: Bacterial recombinant
Tag: No
Stock Buffer: 20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol
Concentration: See tube label
Quality Assurance: >95% by SDS-PAGE

Image



Description: GST is often used as a control in GST pulldown assays.
Storage: Store at -80°C; avoid multiple freeze-thaw cycles
Note: N/A
Literature: N/A