

Agarose 6B

Cat. # P4010-100, P4010-1000

Size:	100 ml (Cat.# P4010-100), 1000 ml (Cat.# P4010-1000)
Bead (Geometry, size):	50 -160 μ m
Cross-Linked:	Yes
Ligand:	N/A
Agarose %:	4% Agarose
Binding Capacity:	N/A
Volume %:	50% (v/v) aqueous suspension containing 20% Ethanol
Application:	Batch, Gravity
Introduction:	Agarose 6B is a bead-formed agarose-based gel filtration matrix. It is also used for coupling affinity ligands to the matrix
Storage Temperature:	2-8 °C

