

PLASTIC 1ML COLUMN

Cat. # P9020

| | |
|-----------------------------|--|
| Introduction: | This column serves as a tool for gravity purification using small quantities of resin. The package contains 5 columns with their top and end caps. The PLASTIC 1ML COLUMN is polypropylene with a total volume of 12 ml and contain a polyethylene frit with a nominal pore size of 20 μm . This column is specially designed for working with bed volumes of 0.5-2.0 ml. |
| Size: | 5 Columns/Pack |
| Column Material: | Polypropylene |
| Frit Material: | Polyethylene |
| Frit pore size: | 20 μm |
| Chemical stability: | Stable in all commonly used reagents |
| Capacity: | 12 ml |
| Caps: | Two caps |
| Storage Temperature: | Room temperature |

