


¹⁵N-labeled Ubiquitin

Cat. # E1130, E1131

Also Known as: ¹⁵N-labeled Ub
NCBI Reference: N/A
MW (no tag): 8.5 kDa
Species: Human
Source: Bacterial recombinant
Tag: No
Stock Buffer: 20 mM Tris, 150 mM NaCl, 2 mM βME, 10% Glycerol
Concentration: See tube label
Quality Assurance: >95% by SDS-PAGE

Image



Description: ¹⁵N-labeled ubiquitin can be used in NMR experiments for determining interactions between ubiquitin and other proteins. It can also be used as a standard for quantification of cellular ubiquitin levels in quantitative mass spectrometric assays.

Storage: Store at -80°C; avoid multiple freeze-thaw cycles

