


6xHis-Sumo 3

Cat. # E3320

Also Known as:	SMT3A; Smt3B; SMT3H1; SUMO-3
NCBI Reference:	NM_006936
MW (no tag):	11.6 kDa
Species:	Human
Source:	Bacterial recombinant
Tag:	6xHis
Stock Buffer:	20 mM Tris, pH7.6 at 4 ⁰ C, 150 mM NaCl, 2 mM βME, 10% Glycerol
Concentration:	See tube label
Quality Assurance:	~95% by SDS-PAGE

Image



Description: SUMO (small Ub-related modifier) is a Ub-like protein. Three types of SUMO are most commonly studied, SUMO 1, SUMO 2, and SUMO 3. SUMO 2 and SUMO 3 are almost identical isoforms and thus share many functions. Like Ub, SUMO can be conjugated to its target proteins as a polymeric chain. However, SUMO 1 forms chains inefficiently as compared to SUMO 2 and SUMO 3. SUMO is conjugated to target proteins by the E1 (SAE1/SAE2), E2 (Ube2I or Ubc9), E3 (RanBP2/Nup358, amongst others). Protein sumoylation is involved in many cellular processes including gene transcription.

Storage: Store at -80°C; avoid multiple freeze-thaw cycles

Note: N/A

Literature:

1. Boddy MN, *et al.* (1996) *Oncogene* 13, 971 – 982.
2. Bayer P, *et al.* (1998) *J Mol Biol* 280, 275 – 286.
3. Melchior F, (2000) *Annu Rev Cell Dev Biol* 16, 591 – 626.
4. Praefcke GJK, *et al.* (2012) *Trends Biochem Sci* 37(1), 23 – 31.
5. Werner A, *et al.* (2012) *Mol Cell* 46(3), 287 – 298.

