


Polyubiquitinated E2-25K

Cat. # G7000

Also Known as: Ubn (K48-linked)-E2-25K
NCBI Reference: N/A
MW (no tag): N/A
Species: Human
Source: Bacterial recombinant
Tag: No
Stock Buffer: 20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol
Concentration: See tube label
Quality Assurance: >90% by SDS-PAGE

Description: Polyubiquitinated E2-25K (Ubn (K48-linked)-E2-25K) is a polyubiquitinated substrate in which various K48-linked polyUb chains were conjugated on the Ub conjugating enzyme E2-25K (GST tagged). It can be used as a substrate to monitor deubiquitination mediated by deubiquitinating enzymes, or to determine binding of K48-linked polyUb chains by proteins that have Ub interacting domains.

Image:



Storage: Store at -80°C ; avoid multiple freeze-thaw cycles

