

PLASTIC 5ML COLUMN

Cat. # P9030

Introduction:	This column serves as a tool for gravity purification using small quantities of resin. The package contains 5 columns with their top and end caps. The PLASTIC 5ML COLUMN is polypropylene with a total volume of 35 ml and contain a polyethylene frit with a nominal pore size of 20 μm . This column is specially designed for working with bed volumes of 2-6 ml.
Size:	5 Columns/Pack
Column Material:	Polypropylene
Frit Material:	Polyethylene
Frit pore size:	20 μm
Chemical stability:	Stable in all commonly used reagents
Capacity:	35 ml
Caps:	Two caps
Storage Temperature:	Room temperature

